Medtronic

PRESSURE DISPLAY BOX

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DRUKWEERGAVE-BOX

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66000EN

Description

The Pressure Display Box is designed to operate for approximately 175 hours using 4 AA (LR6) alkaline batteries, at a temperature range of 15°C to 27°C (60°F to 80°F). The Pressure Display Box has an upper and lower alarm setting. The alarm will be activated when the pressure is equal to or beyond the set value. The alarm settings range from -100 mm Hg to +500 mm Hg.

Indications

This product is intended for use in displaying line pressures during cardiopulmonary bypass surgery.

Contraindications

These devices are not intended for use except as indicated above. Use only Medtronic® DLP Disposable Pressure Display Sets with these devices.

Warnings

The low battery indicator is an important safety feature for a patient monitoring system. A low battery can result in:

- (1) Reduced LCD contrast and therefore a compromised ability to read the LCD display.
- (2) Abnormal behavior in the control of the Pressure Display Box due to low reference levels for logic circuits. Check for this indicator prior to use in any procedure.
- (3) Appearance of all "8"s on the pressure display and all segments on the status display if the battery level is below the required operating voltage on power up. The device will be inoperable until a fresh battery set is installed.

Do not use the device in the presence of a flammable anesthetic mixture with air, or with oxygen or nitrous oxide.

Always use a new, sterile DLP Disposable Pressure Display Set. Never reuse a disposable display set. DLP Disposable Pressure Display Sets should be securely attached before priming. Never attach a display set that is full, or partially full, of fluid to a Pressure Display Box. A full or partially full DLP Disposable Pressure Display Set will reduce the dynamic range and accuracy of the Pressure Display Box.

Never attach a DLP Disposable Pressure Display Set which is under pressure to a Pressure Display Box. Removing the DLP Disposable Pressure Display Set from the Pressure Display Box prior to disconnecting the stopcock may rupture the diaphragm or cause inadvertent fluid loss from the DLP Disposable Pressure Display Set, even at relatively low pressures.

Precautions

Pressure Display Box performance may be adversely affected by electrostatic discharge or electromagnetic interference. Strong fields in all parts of the spectrum may temporarily affect readings or alarm functions. It is recommended that the unit be tested under the conditions and within the environment of normal use to assess these effects.

Pressure Display Boxes contain no user-serviceable parts. Contact Medtronic Cardiac Surgery Customer Service for complete details and return authorization.

Changing the height of the Pressure Display Box relative to the measurement site will require the pressure display to be zeroed.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

NOTE: Proper surgical procedures and techniques are the responsibility of the medical professional. The described procedure is furnished for information purposes only. Each surgeon must evaluate the appropriateness of the procedure based on his or her own medical training and experience, the type of surgical procedure, and the type of monitoring system utilized.

Specifications

Dimensions	12.4 cm x 18.9 cm x 5.1 cm (4.9 in x 7.4 in x 2.0 in)	
Weight	1.1 kg (2.5 lbs)	
Power	4 AA (LR6) (Internally Powered)	
Operation Time	Approximately 175 hours	
Operation Temperature Range	15°C to 27°C (60°F to 80°F)	
Storage Temperature Range	0°C to 50°C (32°F to 122°F)	
Equipment Type	BF	
Range	-100 mm Hg to +500 mm Hg	
Accuracy	-100 mm Hg to 0 mm Hg	(Figure 2)
	1 mm Hg to 300 mm Hg	±5 mm Hg
	301 mm Hg to 500 mm Hg	±10 mm Hg
Moisture Protection	Ordinary	
Mode of Operation	Continuous	
Zero Offset Range	250 mm Hg	


Figure 2.

Adverse Effects

None known.

Directions for Use

- When a Pressure Display Box is first turned ON, both displays will illuminate all segments for 2 seconds as the device tests the display functions. Make a habit of checking the displays each time the device is turned ON. The low battery (LO BAT) indicator on the device should be OFF. Promptly replace batteries whenever the indicator arrow is ON. To install new batteries:
 - Remove battery compartment cover and install 4 AA (LR6) batteries. A good quality alkaline battery type is highly recommended. Rechargeable type batteries (eg, nickel cadmium batteries) are not recommended for use in this unit.
 - Replace cover after installation.
- Securely mount the Pressure Display Box to the structure that will hold it during surgery (eg, bed frame, perfusion cart, etc.). Mount the Pressure Display Box at the same height as the site to be measured (eg, where the catheter tip is placed).
- Attach a DLP Disposable Pressure Display Set, with stopcock open, to the Pressure Display Box. Hold diaphragm housing firmly while rotating collar counterclockwise. Make certain housing is securely tightened.
- With the "FUNCTION" knob in "INSTANTANEOUS" position, zero the display by rotating the "ZERO" knob until the pressure display shows all "0"s. Offsets from at least -100 to +250 are possible.
- With the "FUNCTION" knob in the "ALARM HI" position, the "HIGH" alarm limit can be adjusted using the "ALARM" knob. With the "FUNCTION" knob in the "ALARM LO" position, the "LO" alarm limit can be adjusted using the "ALARM" knob. Both the high alarm set and low alarm set values are displayed on the status display.

NOTE: The high alarm set value cannot be set lower than the lower alarm set value.

6. Activate the desired pressure measurement mode using the "FUNCTION" knob. The mode will also be displayed on the status display.
 - Instantaneous Mode** - Displays pressure waveform. New values are displayed 7 times per second.
 - Systolic Mode** - Displays pressure waveform peaks. Increasing pressures are displayed as soon as a new value of systolic peak is detected.
 - Diastolic Mode** - Displays pressure waveform troughs. Decreasing pressures are displayed as soon as a new value of diastolic peak is detected.
 - Mean Mode** - Displays the mean pressure each time a new value of systolic or diastolic pressure is detected using the following calculation: $[\text{systolic} + (2 \times \text{diastolic})]$ divided by 3.
7. If the pressure being monitored surpasses the pressure measuring range of the Pressure Display Box, the status display will flash the upper or lower pressure limit values (+500 mm Hg or -100 mm Hg, respectively) and the alarm warnings will activate.
8. The alarm will be activated when the pressure value equals or exceeds the alarm set values. The alarm warns the user with a flashing red light, a flashing high alarm set value or the low alarm set value on the status display and an audible alarm. Setting the alarm toggle switch into the "OFF" position will mute the audible alarm. The alarm functions will be active until the alarm condition is corrected.
9. Following surgery, remove the DLP Disposable Pressure Display Set by disconnecting the distal stopcock first and isolating the Pressure Display Box from the patient. Detach the pressure dome from the Pressure Display Box by rotating the threaded collar clockwise.
10. If the Pressure Display Box will not be used for 30 days or more, remove the batteries from the unit.
11. Clean and disinfect the Pressure Display Box with a soft cloth and cleaning/disinfecting detergent.

EQUIPMENT LIMITED WARRANTY¹

The following limited warranty applies to United States customers only:

- A. This LIMITED WARRANTY provides the following assurance to the purchaser of the Pressure Display Box (model 66000EN), hereafter referred to as the "Equipment":
- (1) Should the Equipment fail to function within normal tolerances due to a defect in materials or workmanship within a period of one (1) year, commencing with the delivery of the Equipment to the purchaser, Medtronic will at its option: (a) replace any defective part or parts of the Equipment; (b) issue a credit to the purchaser equal to the Purchase Price, as defined in Subsection A(2), against the purchase of the replacement Equipment; or (c) provide a functionally comparable replacement Equipment at no charge.
 - (2) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement, Equipment.
- B. To qualify for this replacement or credit set forth in Section A, the following conditions must be met:
- (1) The Equipment must be returned to Medtronic within thirty (30) days after discovery of the defect.
 - (2) The Equipment must not have been repaired or altered outside of Medtronic's factory in any way which, in the judgment of Medtronic, affects its stability and reliability. The Equipment must not have been subjected to misuse, abuse or accident.
- C. This LIMITED WARRANTY is limited to its express terms. In particular:
- (1) Except as expressly provided by this LIMITED WARRANTY, MEDTRONIC IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE OR MALFUNCTION OF THE EQUIPMENT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.
 - (2) This LIMITED WARRANTY is made only to the purchaser of the Equipment. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOM OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PATIENT SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.
 - (3) The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid. This LIMITED WARRANTY gives the purchaser specific legal rights. The purchaser may also have other rights which vary from state to state.
 - (4) No person has any authority to bind Medtronic to any representation, condition or warranty except this LIMITED WARRANTY.

¹. This limited warranty is provided by Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432. Areas outside the United States should contact their local Medtronic representative for exact terms of the Limited Warranty.

EQUIPMENT LIMITED WARRANTY¹

The following limited warranty applies to customers outside the United States:

- A. This LIMITED WARRANTY provides the following assurance to the purchaser of the Pressure Display Box (model 66000EN), hereafter referred to as the "Equipment", that should the Equipment fail to function within normal tolerances due to a defect in materials or workmanship within a period of one (1) year, commencing with the delivery of the Equipment to the purchaser, Medtronic will at its option: (a) replace any defective part or parts of the Equipment; (b) issue a credit equal to the original Equipment purchase price (but not to exceed the value of the replacement Equipment), against the purchase of replacement Equipment; or (c) provide functionally comparable replacement Equipment at no charge.
- B. To qualify for this replacement or credit, the following conditions must be met:
 - (1) The Equipment must be returned to Medtronic within thirty (30) days after discovery of the defect.
 - (2) The Equipment must not have been repaired or altered by someone other than Medtronic in any way which, in the judgment of Medtronic, affects its stability and reliability.
 - (3) The Equipment must not have been subjected to misuse, abuse or accident.
- C. This LIMITED WARRANTY is limited to its express terms. In particular, Medtronic is not responsible for any incidental or consequential damages based on any use, defect or failure of the Equipment, whether the claim is based on warranty, contract, tort or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

¹ This limited warranty is provided by Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432. Areas outside the United States should contact their local Medtronic representative for exact terms of the Limited Warranty.

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Medtronic of Canada Ltd.